

TECHNICAL SPECIFICATION (General)

Total No. of Outputs	Two
Number of Channels	Four Channel 1 & 2 for IFT/Galvanic/Russian current mode Channel 3 & 4 for Faradic current mode
Cyclic ON	1-30 sec. selectable
Cyclic OFF	0-30 sec. selectable
Treatment Time	1-60 minutes selectable
Power	110VAC, 60 Hz or 220VAC, 50 Hz
Operating Temperature	32 F - 112F
Stimulation types	Interferential/Russian/Galvanic & Faradic Stimulation
Weight	4.22 lb
Dimensions	12.2" (l) x 3.5" (h) x 10" (b)
Treatment timer indicator	The Digital window indicates the remaining treatment time in every one minute, for each output.
Ramp Rise Time	Two Seconds
Ramp Fall Time	Two Seconds

Note : Cyclic ON & Cyclic Off not applicable in Vector. When Cyclic On & Cyclic Off key are pressed the digital display goes "---" & buzzer beeps.

In dual mode when pulse rate key is pressed the digital display goes "---" & buzzer beeps.


 johari[®]
digital