

Digi Tens - J249 D

johari[®]
digital

TECHNICAL SPECIFICATION

No. of Output	:	Two
Stimulation	:	Symmetrical Monophasic square wave
Mode	:	Selectable burst, conventional and modulation
Treatment Time	:	Selectable 1-60 minute and continuous
Power	:	9V Alkaline battery & 9V Adaptor
O/P Voltage	:	40 Volt @ 500 Ohm
O/P Current	:	100 m Amp @ 500 Ohm load
Pulse Width	:	50 to 260 μ s in step of 10 μ s
Pulse Rate	:	002 Hz to 150 Hz, in step of 1 Hz upto 20 Hz, in step of 5Hz from 20 Hz to 100 Hz In step of 10 Hz from 100 Hz upto 150 Hz
Display System	:	Customized LCD

johari[®]
digital