

TECHNICAL SPECIFICATION

Channel	Two
Carrier Frequency	4000 Hz ($\pm 10\%$) Fixed
Difference Frequency	1 ($\pm 3\text{Hz}$) to 250 ($\pm 10\%$) Adjustable.
Output Voltage	27V \pm 2.5 Volts (500 Ohms load)
Power Source	Battery 9V Alkaline and Adaptor
Stimulation types	Interferential Stimulation
Waveform	Symmetric Biphasic Square wave
Four Frequency Shift	Continuous, 1\1 Abrupt, 8\8 Abrupt, 10\10 Ramped Options