

TECHNICAL SPECIFICATION

Number of Output	Four
Number of Channels	Four
Stimulation Type	EMS
Waveform	Biphasic
Output Voltage	50 V PP @ 500 Ohm Load
Output Current	100mA @ 500 Ohm Load
Pulse Frequency	Selectable Burst 200, 600 and 800 Hz
Pulse width	0.5 milliseconds
Cyclic Off	1 - 60 Sec. (Selectable)
Cyclic On	0 - 60 Sec. (Fixed)
Timer	15, 30, 45, 60 min. (Selectable)
Power Source:	90-250V AC
Temperature	32 F 112 F
Weight	5.1 lb
Dimensions	12.2 "(1 x 3.5" (h) x 10" (b)