

TECHNICAL SPECIFICATION

Power supply : 220V AC, 50Hz or 110VAC, 60 Hz
 Wave Shape : biphasic Square Wave
 No. of Channel : 4
 Modes : 7
 Output Voltage : 50 V @ 560 Ohm

MODES	PULSE WIDTH	FREQUENCY
Lifting	300 Micro Sec.	200 Hz.
Firming	300 Micro Sec.	100 Hz.
Toning	300 Micro Sec. In 6 steps.	100 Hz.
Reshaping	300 Micro Sec.) In 5 steps.	130 Hz.
Massaging	300 Micro Seconds	16 Hz.
Pain Relief	200 Micro Seconds	50 Hz.
Auto All	Mode changes after Every 1 minutes.	

