

TECHNICAL SPECIFICATION

Total No. of Outputs	Four
Total No. of Channels	Eight through Four Outputs
Output Power	0 - 60V PP @ 500 ohms Load
Cyclic ON	1-60 Sec. Selectable (not applicable in IFT Mode)
Cyclic OFF	0-60 Sec. Selectable (not applicable in IFT Mode)
Timer	1-60 minutes selectable
Power	220VAC / 50Hz or 110VAC / 60 Hz
Temperature	32F-112F
Stimulation types	Interferential and Russian
Dimensions	45 cm(l) x 16 cm (h) x 34 cm (b)

